

<i>Analytická osobní váha • NÁVOD K OBSLUZE</i>	CZ	4-12
<i>Analytická osobná váha • NÁVOD NA OBSLUHU</i>	SK	13-21
<i>Analytic bathroom scale • INSTRUCTIONS FOR USE</i>	EN	22-29
<i>Analitikus személyi mérleg • KEZELÉSI ÚTMUTATÓ</i>	HU	30-37
<i>Osobowa waga analityczna • INSTRUKCJA OBSŁUGI</i>	PL	38-45
<i>Analytische Personenwaage • BEDIENUNGSANLEITUNG</i>	DE	46-55

TRACEY

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HOUSEHOLD USE ONLY. DO NOT IMMERSE IN WATER OR OTHER LIQUIDS. OPEN.

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TO AVOID DANGER OF SUFFOCATION, KEEP THIS PLASTIC BAG AWAY FROM BABIES AND CHILDREN. DO NOT USE THIS BAG IN CRIBS, BEDS, ARRIAGES OR PLAYPENS. THIS BAG IS NOT A TOY.

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TRACEY

INSTRUCTIONS FOR USE

Dear customer, thank you for purchasing our product. Please read the operating instructions carefully before putting the appliance into operation and keep these instructions including the receipt and, if possible, the box with the internal packing.

I. SAFETY WARNING

- Consider the instructions for use as a part of the appliance and pass them on to any other user of the appliance.
- **The appliance is intended for use in households only and for similar use! It is not intended for use in health-care facilities and for commercial use!**
- This appliance can be used by children aged from 8 years and above and persons with reduced physical, sensory or mental capabilities or lack of experience and knowledge if they have been given supervision or instruction concerning use of the appliance in a safe way and understand the hazards involved. Children shall not play with the appliance. Cleaning and user maintenance shall not be made by children without supervision.
- Never use the appliance if it malfunctions, if it fell to the ground and got damaged or if it fell into water. In these cases bring the appliance to a specialised service to have the safety and functions of the device checked.
- **We recommend: Do not step onto the scale with wet feet or on wet surface of the scale, there is a risk of slipping.**
- **Never immerse the scale into water (even their parts)!**
- Keep in mind that the data read from this scale are approximate and inaccurate in comparison to results of actual medicinal analyses. The exact determination of biometrical characteristics can be done only by a doctor via medical methods.
- Do not use during pregnancy.
- Do not stand on the edge of the scales!
- Do not use the fat analysis feature on this scale if you are using one of the following devices: pace maker, electrocardiograph and/or other medical device implanted in the body or device used as support.
- If the scale was stored in lower temperatures, let it adjust to the room temperature.
- Use the scale in a horizontal position only and at places where they cannot knock over, in sufficient distance from heat sources (e.g. **stove, fireplace, heat radiators**), sources of increased air humidity (e.g. **sauna, bathroom, swimming pool**) and equipment with strong electromagnetic field (e.g. **microwave oven, radio**).
- Protect the personal scale from dust, chemicals, excessive humidity, toiletries, liquid cosmetic preparations, intensive temperature changes and do not expose them to direct sunlight.
- Weigh yourself without clothing and footwear, before meals and always at the same time of day.
- Do not place the scale on uneven surface or on surface covered with carpet. Unstable or soft surface under the scale can have negative effect on precision of weighing.

- In handling the scale, make sure injury of persons or damage of the appliance is prevented.
- Do not overload the appliance weight greater than the maximum weight capacity! This prevents damage to the appliance.
- Handle the scale care. Do not throw it or jump on it. The scale is durable, but improper handling can destroy its electronic sensors.
- Do not put any objects on the surface of the scale during storing, they could be damaged.
- Do not use rechargeable accumulators.
- If the battery is leaking, remove it immediately; it can damage the scale.
- Dispose of a discharged battery in a suitable way (see par. **IX. ENVIRONMENTAL PROTECTION**).
- Keep batteries and scale out of reach of children and incapacitated people.
A person, who would swallow the batteries, must seek medical help immediately.
- Never use the appliance for any other purpose than for that described in these instructions for use!
- **WARNING:** There is a risk of injury in case of incorrect use (not in accordance with the manual) of the appliance.
- The producer is not responsible for damage caused by improper use of the appliance (e.g. **breaking the board, etc.**) and its does not cover the appliance in the case of failure to comply with the safety warnings above.

Important notes

- This scale requires you to weigh yourself barefoot for the reason of correct conductivity during the body fat analysis.
- Do not place the scale on an uneven surface or a surface covered with a carpet.
- Excessive drinking, eating, exercising, health conditions, used medicaments, women's menstruation cycle etc. influence the results of the analysis.
- Here must not be a contact between feet, calves or thighs. Otherwise the measurement could be inaccurate.
- Simply stepping on the weighing surface is not enough. Grasp the sensor with the whole palm of your hand as shown in the figure below. Grasp the left arm with the left palm and the right arm with the right palm.
- Due to the inherent differences between people, the instructions are for information purposes only.

II. GENERAL INFORMATION

You have purchased an advanced device that is a multifunctional modern assistant for monitoring and analyzing various body characteristics.

Body Fat Ratio and Body water Ratio Chart

Age	Female		Male		Result
	Body fat (%)	Body water (%)	Body fat (%)	Body water (%)	
< = 30	4.0 - 20.5	66.0 - 54.7	4.0 - 15.5	66.0 - 58.1	Underfat
	20.6 - 25.0	54.6 - 51.6	15.6 - 20.0	58.0 - 55.0	Healthy
	25.1 - 30.5	51,5 - 47.8	20.1 - 24.5	54.9 - 51.9	Overfat
	30,6 - 80,0	47,7 - 13,8	24,6 - 80,0	51,8 - 13,8	Obese

< 30	4.0 - 25.0	66.0 - 51.6	4.0 - 19.5	66.0 - 55.4	Underfat
	25.1 - 30.0	51.5 - 48.1	19.6 - 24.0	55.3 - 52.3	Healthy
	30.1 - 35.0	48.0 - 44.7	24.1 - 28.5	52.2 - 49.2	Overfat
	30,6 - 80,0	47,7 - 13,8	24,6 - 80,0	51,8 - 13,8	Obese

Bone Weight Chart

Weight (kg)	Female			Male		
	< 45	45 - 60	> 60	< 60	60 - 75	> 75
Ref. bone weight (kg)	> 1,8	> 2,2	> 2,5	> 2,5	> 2,9	> 3,2

Musculature content Chart

Age	Female %			Male %		
	Low	Normal	High	Low	Normal	High
<30	<34	34 - 43	>43	<42	42 - 57	>57
31 - 50	<31	31 - 38	>38	<40	40 - 52	>52
51>	<27	27 - 34	>34	<37	37 - 48	>48

KCAL Chart

Age	Female			Male		
	Weight (kg)	Basal metabolic value (kcal / kg / day)	KCAL	Weight (kg)	Basal metabolic value (kcal / kg / day)	KCAL
15 - 17	51,4	25,3	1300	59,8	27	1610
18 - 29	51,2	23,6	1210	64,7	24	1550
30 - 49	54,2	21,7	1170	67	22,3	1170
50 - 69	53,8	20,7	1110	62,5	21,5	1110
>= 70	48,7	20,7	1010	56,7	21,5	1010

BMI Chart

Symbol				
Status	Underfat	Healthy	Overfat	Obese
BMI	< 18,5	18,5 - 27	30 - 39	> 40

Due to the natural differences between people this chart should be used for reference purposes only. For the reason of differences between people, this table is only informative. Always consult your attending physician if any measuring results are abnormal. The results of measuring bones' weight and muscular proportions are simply tentative. Healthy man has approximately 40 % of muscular proportion and healthy woman has 30 % according to an average ratio. Women's weight of bones is among 12 % and men's weight of bones is among 15 % but it depends on many relevant factors.






III. FEATURES OF THE DEVICE (FIG. 1)

A – Unit of scale

- A1 – display
- A2 – glass surface
- A3 – measuring electrodes
- A4 – SET button  – power on, setting and confirmation

- A5 – button  and  = selection buttons
- A6 – UNIT button – change of unit mode (kg, st:lb, lb, ft, cm)
- A7 – battery cover

LCD SEGMENTUM / ICONS ILLUSTRATIONS:

	Kg / lb / st:lb weight measurement units
	P0 ...P9 number of memory
	 man / woman
	MUS musculature contents (%)
	BONE weight contents (%)
	cm (ft: in) height measurement units
	Age age
	FAT body fat contents (%)
	Underfat / Healthy / Overfat / Obese
	TBW body water contents (%)
	BMI body mass index
	KCAL minimum energy consumption

IV. PREPARATION FOR USE

Remove all the packing material and take the scale out. Remove all possible adhesion foils, stick-on labels or paper from the scale. Open the battery cover **A7** on the bottom of the scale, insert the battery **with the correct polarity** and close the cover in the opposite direction. Place the scale on flat, stable and dry surface. Switch on the scale by tapping on the board of the scale with your foot and check the displayed weight unit (**kg**). If you wish to set unit to (kg, st:lb, ft, cm) press the UNIT button on the bottom of the scale (Fig. 3)

Weight and height units

- kg = kilogram
- lb = pound (British weight unit)
- st:lb = stone:pound (American/British unit of weight)
- ft, cm = feet, centimetre

- 1 kg = 2.2 lb
- 1 lb = 0.454 kg
- 1 st = 6.356 kg
- 1 FT = 30.48 cm

V. INSTRUCTIONS FOR USE

Your scale is fine electronic equipment. One of the things that are most difficult to weigh is human body as it is always in motion.

If you want to obtain reliable results, always try to step with your feet at the same place on the scale and stand still. If possible, weigh yourself at the same time of a day (preferably in the morning) to obtain comparable results. Remember that excessive clothes can increase your weight as well as eating and drinking before weighing. Short-term weight fluctuation is generally attributed to loss of liquids.

A) Regular weighing mode

Turn the scale on by stepping on the weighing surface. After a moment, the measured weight flashes and the value remains displayed for several seconds. **Then the scale are switched off automatically.**

Note

- The scale is switched on automatically if you do not step on it within approx. 10 seconds after „0.0“ is displayed.
- We recommend using the scale from about 6 kg of the subject.

B) Weighing and total body analysis mode

Setting (change) of personal data for the body fat and water analysis

Before the first analysis, it is necessary to put in personal data (sex, height and age) which will be saved to the memory. It is possible to save data of 10 persons in the scale's memory. In order to ensure accurate results users should enter the required personal parameters; age, sex and height. These scale have the facility to store the personal parameters of 10 people.

The input data can be set and gradually changed in the following range:

Memory	Gender	Height	Age
0 – 9	male / female	100 – 240 cm / 3'3" – 7'10"	10 – 80

SETTING AND SAVING PARAMETERS

Press **SET** to turn ON.

Press  or  to switch between Body Parameter P0 - P9. Press **SET** to confirm.

Press  or  to select male or female. Press **SET** to confirm.

Press  or  to indicate age. Press **SET** to confirm.

Press  or  to indicate height. Press **SET** to confirm.

Now all the personal data is saved (changed). The values are flashed and recorded.

Minimum energy consumption CAL

It represents optimum amount of energy consumption necessary for maintaining the basic vital functions and balanced metabolism. The data on the display show the optimum daily consumption of energy that you give out at your weight, height, age and gender. You should also take in this amount of energy to keep your metabolism balanced. If you intend to reduce your weight, your daily energy intake should be a bit lower. However, remember to keep carbohydrates, fats and proteins balanced in your diet. 239 cal = 1 kJ

Important:

During body fat testing, please remove shoes and socks, etc., and clear your feet, in order to get the most accurate result. To ensure accuracy, please test at the same time of the day.

Test result may be inaccurate when the users are:

- after intense sport activity, overfeeding, or excessive dehydration,
- children under 10, or adults over 80,
- body building athletes, or other occupational athletes,
- pregnant women; Dialysis or edema patient.

MEASURING

- 1) Tap on the scales platform to switch on the product
- 2) Press **SET** P1 will flash. Press  or  to switch between Body Parameter P0 - P9.
- 3) Step onto the platform when the LCD display reads 0.0.
Important: For correct measurements and scale functions step on the scale barefoot.
- 4) First the weight is recorded, then „  “ is displayed, indicating that the scales are analyzing. Please do not alight.
- 5) Once finished, the weight, body fat, body muscle, bone weight, BMI and KCAL will be displayed alternately for 2 times before the scales automatically turn off.
- 6) After this the scale is going turn off.

Compare the body fat reading against the <Body Fat Radio Chart> to confirm the body status.

Replacing batteries

Open the battery cover **A7** on the bottom of the scale (Fig. 2).. Insert a battery with the right polarity and close the cover. Use AAA batteries (3 pcs) for this appliance. If you are not using the appliance for some time, take the batteries out!

Notes

- Do not combine various types of batteries. Do not even combine used and new batteries together.
- Batteries not intended for this purpose must not be charged.
- Remove the batteries if you will not be using the appliance for a long time.
- The supply terminals must not be short-circuited.

VI. TROUBLESHOOTING OF POSSIBLE PROBLEMS

If you have problems with weight:

- Check whether the battery is inserted properly.
- Check whether you selected the proper weight unit.
- Check whether the scale is on a flat floor and they do not touch the wall.
- If nothing appears on the scale when you want to switch them on by tapping on them or if „**Lo**“ appears only, replace the batteries.
- If „**Err**“ appears on the display, the scale was overloaded.
- If „**Err2**“ appears on the display, the scale indicates an error in measuring the fat percentage.
- If „**FAL**“ / „**FAH**“ appears on the display, it shows an error in the body fat percentage. Fat percentage exceeds the maximum capacity is over 50 %, fat percentage is under 5 %.
- If unusual signs are displayed, remove the insert the battery again.
- If signs appear, such as if the display is flashing or if there is a faulty sign, try to move the scale away from the influence of a source of interference or ensure switching off the source for the time of using the scale.

- Abnormal measurement data. Please do not wear shoes and socks for measurement, and remove metal accessories and other metal jewel from your body.

VII. STORING

After using the scale, put it into the horizontal position (not vertically), out of reach of children and incapacitated people.

VIII. MAINTENANCE


Do not use coarse and aggressive cleaning agents (e.g. **sharp objects, scrapers, diluting agents or other solvents**). Clean the board by wiping with a soft wet cloth. It could be slippery when the unit platform surface is wet, or semi-damp. Ensure the weighing platform surface is dry before using for safety reason. Make sure that water does not get inside the appliance.

More extensive maintenance or maintenance requiring intervention in the inner parts of the appliance must be carried out by professional service!

Manufacturer contact: info@eta.cz

Failure to comply with the instructions of the manufacturer will lead to expiration of the right to repair!

The information regarding warranty and product repairs you can find on the website www.eta.cz/servis-eu.

IX. ENVIRONMENTAL PROTECTION

If the dimensions allow, there are labels of materials used for production of packing, components and accessories as well as their recycling on all parts. The symbols specified on the product or in the accompanying documentation mean that the used electrical or electronic products must not be disposed of together with municipal waste. For proper disposal of the product, hand it over at special collection places where it will be taken over free of charge. Suitable disposal of the product can help to maintain valuable natural resources and to prevent possible negative impacts on the environment and human health, which could be the possible consequences of improper waste disposal. Fines can be imposed for improper disposal of this type of waste in agreement with the national regulations. Remove a discharged battery and dispose of it in a suitable and safe way through the special collection networks. Never dispose of the battery by burning!

X. TECHNICAL DATA

Weighing range max.:	180 kg / 396 lb / 28 st
Discrimination:	100 g
Net weight of device:	1.55 kg
Size of the product (mm):	300 x 300 x 23 mm
Protection class of the appliance:	III.

The manufacturer reserves the right to change the technical specifications for the respective model.

WARNINGS AND SYMBOLS USED ON THE APPLIANCE, PACKAGING OR IN THE
INSTRUCTIONS MANUAL:

HOUSEHOLD USE ONLY. DO NOT IMMERSE IN WATER OR OTHER LIQUIDS.

TO AVOID DANGER OF SUFFOCATION, KEEP THIS PLASTIC BAG AWAY FROM BABIES AND CHILDREN. DO NOT USE THIS BAG IN CRIBS, BEDS, CARRIAGES OR PLAYPENS. THIS BAG IS NOT A TOY.

WARNING

It is not allowed to dispose it together with communal waste

Please read the instructions manual

Do not step onto the scale with wet feet or on wet surface of the scale, there is a risk of slipping.

AUTO STEP ON – just step on the scale and stand still, the scale will automatically turn ON and start taking measurements.

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Service: Omega electric GmbH, Servicezentrale, Industriering 2, D- 04626 Schmölln.